
caring mill™

TENS Pain Management Solution

Model 25517 (VH 56-032)

INSTRUCTION MANUAL ENGLISH & ESPAÑOL

Intention for Use:


The TENS Pain Management Solution is intended for temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities.

PLEASE READ THESE INSTRUCTIONS
COMPLETELY BEFORE USING YOUR TENS DEVICE.


STOP! PLEASE ENSURE YOU HAVE ALL OF THE FOLLOWING COMPONENTS BEFORE USING YOUR TENS DEVICE.

Control Unit

1 Set of Gel
Electrode Pads

Storage Case

One CR2032
Battery

Instruction Manual
Replacement Electrode Order Form

**IF YOU ARE MISSING ANY PARTS, INCLUDING INSERTS OR INSTRUCTION MANUALS, DO NOT RETURN TO PLACE OF PURCHASE.
CONTACT CUSTOMER CARE AT 866-326-1313.**

INDEX

Care and Safety Information.....	4-11
Product Features and Symbols.....	12-13
Battery Installation / Replacement	14
Preparing Your Device	15
Placement & Treatment.....	16-17
Instructions for Use.....	18-21
Cleaning and Storage.....	22
Troubleshooting	23
Device & Label Symbols	24
Electromagnetic Compatibility.....	25-28
FCC Statement.....	29-30
Product Specifications	31
Warranty Information	32-33
Español.....	35-68

Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

Caring Mill™ 2017
240 W. 37th St.
6th Floor
New York, NY 10018
1-888-372-1450
Made in China

#93-1387 08/17

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

NOTE: Read all instructions and inserts included with this device carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in injury or equipment damage.

Improper care or use of your device may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the device's efficacy and long life.

CONTRAINDICATIONS :

Do not use this device if you have a cardiac pacemaker, implanted defibrillator or other implanted metal or electronic device.

Do not use this device to treat undiagnosed chronic pain.

Do not use this device if you have suspected or diagnosed heart problem.

Do not use this device if you have suspected or diagnosed epilepsy.

Do not use this device if you have a tendency to bleed internally following an injury.

Do not use this device if you recently had surgery, or have ever had surgery on your back.

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

Do not use this device if areas of skin lack normal sensations, such as skin that tingles or is numb.

Do not use this device during menstruation or during pregnancy.

GENERAL WARNINGS PRIOR TO USE

- If you are under the care of a Physician, consult with your Physician before using this device
- Long-term effects of this device are not known
- Do not place the device on or close to your heart
- Take care when applying device close to the neck; do not place this device on the front or sides of the neck. Severe spasm of the muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing. Stimulation over the neck could also affect hearing or blood pressure
- Do not apply stimulation across the chest because the introduction of electrical current into the chest may cause rhythm disturbances to the heart
- Do not place the device on or around your head. The effects of stimulation of the brain are unknown

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- Do not use this device while sleeping
- Do not use if you feel numbness
- Do not use this device in or close to water
- Use the only on normal, healthy, clean and dry skin. Do not use on or close to open wounds or rashes, or over swollen, red, infected or inflamed skin
- Do not apply stimulation over, or in proximity to, cancerous lesions
- Do not use the device on children or incapacitated persons
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals

GENERAL CAUTIONS PRIOR TO USE

- Read this User Manual before using this device for the first time
- Keep this manual available whenever you use this device
- This device is intended for individual person use only
- This device is not effective for pain associated with Central Pain Syndromes, such as headaches

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- This device is for pain caused by muscle soreness, and should be placed only around muscles where pain originates
- The pain may indicate that you have some other health problem. You should know the reason and source of your pain before using these devices. Do not solely rely on the treatment of this device for pain
- The safety of using this device during pregnancy or birth has not been established
- The effectiveness of this device depends greatly on a person's individual physical condition. It may not always be effective for every user
- If you have had medical or physical treatment for your muscle pain, consult with your treatment provider before using this device. You should contact your physician prior to using this device following recent surgical procedures. Stimulation may disrupt the healing process

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

OPERATING WARNINGS & CAUTIONS

- CONSULT WITH YOUR PHYSICIAN IF YOU HAVE ANY OF THE FOLLOWING:
- If you have suspected or diagnosed heart problem
- If you have suspected or diagnosed epilepsy
- If you have a tendency to bleed internally following an injury
- If you recently had surgery, or have ever had surgery on your back
- If areas of skin lack normal sensations, such as skin that tingles or is numb
- During menstruation or during pregnancy
- The unit is intended for the temporary relief of pain caused by muscle soreness. The source or cause of the pain should be addressed by a healthcare professional. Pain could indicate an injury or condition that requires treatment
- Some people may feel skin irritation or experience a very sensitive feeling in the skin due to electrical stimulation. If this occurs, stop using these devices and consult your physician
- If skin under the pad feels irritated after using the stimulator for a long period of time, use the stimulator for a shorter period of time

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- Minor redness at stimulation placement is a normal skin reaction. It is not considered a skin irritation, and it will normally disappear within 30 minutes after the electrodes are removed. If the redness does not disappear after 30 minutes from the removal of electrodes, do not use the stimulator again until after the excessive redness has disappeared
- Turn off your device if the stimulation feels unpleasant or does not provide pain relief
- Keep this device out of the reach of children.
- Use this device only with the pads and accessories recommended by the manufacturer
- Do not use this device when driving, operating machinery or when swimming
- Before removing the Electrode Pad, be sure to turn it OFF, avoiding unpleasant stimulation
- TENS is not a substitute for pain medications and other pain management therapies
- TENS devices have no curative value

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS

- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism
- Do not use the unit while sleeping
- After each use, clean any adhesive residue left on the skin with soap and water
- The effectiveness of stimulus therapy varies from person to person. It may not be an effective treatment for all users
- The long-term effects of stimulus therapy are unknown

STORAGE WARNINGS & CAUTIONS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: 14°F – 140°F (10°C – 60°C); humidity: 30-95% RH
- Keep the unit out of reach of small children
- Remove the batteries if the unit will not be used for an extended period of time

IMPORTANT SAFEGUARDS WARNINGS & CAUTIONS


- Do not store the device in direct sunlight, dusty or humid environments, or extreme temperatures. Exact figures for appropriate use and storage can be found in the Product Specifications section of this manual

CLEANING WARNINGS & CAUTIONS

- Never immerse the unit in water to clean as it may damage the unit
- Follow the 'Cleaning and Maintenance' portion of this manual for instruction on how to clean and care for your device

PRODUCT FEATURES

Control Unit

TENS
Control Unit

Increase
Intensity

+


Decrease
Intensity

-

Power
Button

PRODUCT FEATURES

Gel Electrode Pads




BATTERY INSTALLATION

This unit comes with 1 C2032 battery pre-installed. To begin using, carefully pull the protective tab from the battery compartment; the device is now ready for use.

To replace the batteries:

1. Use the included 'Battery Compartment Opener' or a coin to open the battery cover.
2. Follow the indicators on the battery cap—left/open; right/lock—to open the battery cap. If necessary, use a paper clip to pop off the cap.
3. Remove the expired battery; you may need a pin or paper clip to pop out the battery.
4. Replace C2032 battery into the battery compartment, with the plus (+) side facing up. Always use a new battery.
5. Replace the battery cover.
6. Dispose of batteries according to local disposal and recycling regulations.




It is recommended to remove the battery if the unit will not be used for an extended period of time.

PREPARING YOUR DEVICE

It is necessary to prepare and assemble your device prior to treatment.

NOTE: Save the large plastic film for storage after treatment to protect the gel pads from dirt and debris.

1. Peel off the one side of the blue plastic from the gel pads.
2. The gel pads are applied to the back of the TENS Control Unit along the curved electrode ridges.
3. While the device is OFF (the GREEN light is NOT illuminated), carefully align the gel pads on the back side as shown and press firmly to ensure the pads are completely adhered to the electrode ridges.
4. Carefully remove the plastic film from the electrode pads.
5. You are ready to apply the pad to the affected treatment area.

Do not turn the power ON until the device is in place on the user's skin as indicated in the 'Instructions for Use' section of this manual.

PLACEMENT & TREATMENT


Read all warnings and cautions completely before use.

The electrode pads are applied directly to the skin. Clothing may be worn over the device; clothing should not constrict the device and should allow access to the control unit buttons.

Targeted treatment areas may include:

- shoulder
- elbow
- thigh
- knee
- calf
- ankle
- back of neck
- upper back
- lower back

Never place device on over the heart, on the front or sides of neck or the face.


PLACEMENT & TREATMENT

The TENS device may be placed directly on the treatment site where you are experiencing pain.

An alternative therapy is to apply the TENS device on alternating sides of the location of the pain.


For example, if the treatment site is the elbow, place the TENS unit above the pain site for the 20 minute cycle, and then place below the pain site for 20 minutes, surrounding the location.

INSTRUCTIONS FOR USE

Read all warnings and cautions completely before use.
Treatment time is 20-minutes. Discontinue use immediately if treatment is uncomfortable or painful.

1. Treatment area should be clean and dry before applying the control unit. Note that excessive skin lotion or soap residue may affect the adhesiveness of the pad. Only use on HEALTHY skin areas with no cuts, sores or wounds.
2. Apply the control unit and pad to the desired treatment area.
3. Ensure the electrode pads are smoothly applied and the gel adhesive is holding the pad firmly in place on the skin.
4. Press the POWER symbol button for 3 seconds. A GREEN led light will illuminate to indicate the power is on.
5. The TENS unit will begin to operate at the lowest level setting. The default lowest level mode is a rotation of 4 pulse types at level 1 — scraping, massage, acupuncture and tapping.
6. Press and release the ‘+’ button to find your desired therapy sensation level. This device offers 15 intensity levels. Press and release the ‘-’ button to reduce the intensity level.

INSTRUCTIONS FOR USE

As you press and release the button, the green light will illuminate once for each intensity level adjustment.

7. To change the mode, while the unit is on/operating, press and hold the '+' for 3 seconds; the mode will switch the next in the sequence of pulse types — scraping, massage, acupuncture, tapping or combination (all 4 types in sequence/default mode). As you press and hold the button, the green light will illuminate two times to indicate the change in mode type.
8. The sensation should be strong but comfortable; press the + or – button at any time during the treatment to change the intensity as desired.
9. The device will pulse for several seconds, pause and resume repeatedly until the treatment is finished.
10. The device will automatically shut off after 20 minutes of use. To discontinue use prior to the end of treatment time, press and hold the POWER button. The green light will flash three times and the device will turn off. Do NOT attempt to remove the device during operation.

Until you are familiar with use of the unit, begin with the lower intensity levels and gradually progress to more advanced levels as you feel comfortable.

INSTRUCTIONS FOR USE

11. Be sure the device is OFF (no sensation is felt) before removing the device from your skin. Take care during removal to best protect the adhesive and preserve the life of the electrode pad.
12. Place the plastic cover over the gel pads to prevent them from drying out. Store the device in the included storage case.
13. To conserve the battery and prevent shocks, the device will automatically power off when not in use. The green light will illuminate for three seconds, then flash three times before shutting off.
14. It is recommended to remove the batteries if the unit will not be used for an extended period of time

INSTRUCTIONS FOR USE

Note:

You may wish to clean the treatment site after use to remove any adhesive residue from the electrode pad.

Gel electrode pads

- The electrode pads can be reused as long as they retain their stickiness. Cover the pads with the included film pieces when not in use to prevent them from drying out
- Properly cared for, the electrode pads may last up to 50 uses; the life of the electrode pad is largely dependent upon proper care and the condition of the skin as applied for use
- To help prolong the life of the electrode pad, between uses, moisten the gel side of the pad with a drop of water and allow to dry; replace the film for storage
- If you notice that the pads have become dry, no longer adhere well to the skin, or begin to irritate the skin, they are in need of replacement.

Replacement pads are available through **FSAstore.com**.

Please use model number 25518.

CLEANING AND STORAGE

Cleaning

- Remove the battery before cleaning the device
- The Control Unit can be cleaned with a soft, dry cloth
- Never use cleaning agents or excessive water to clean any part of the unit. Do not submerge any part of the unit or let it get wet
- Never disassemble the unit or attempt to repair it (other than replacing the adhesive pads and batteries)

Storage

- The gel pads must remain on the control unit between uses; place the film cover over the electrode pads to prevent drying out or getting dirt and debris on them
- Protect the unit from mechanical shock or heavy impact
- Avoid exposure to extreme temperatures
- Remove the batteries during extended storage

TROUBLESHOOTING

Problem	Solution
LED does not illuminate / Control unit will not turn on	Check the polarity of the battery. Replace the battery.
LED illumination is weak	Battery power is low, replace the battery.
Intensity output is low / no stimulation is occurring	Electrode pad may be dirty, clean and attempt again. Confirm that the control unit is firmly attached to the electrode pads and attempt again.
Device turns off prematurely / before the 20-minute time is up	Battery power is low, replace the battery. Electrode pad may be dirty, clean and attempt again Device may not be applied properly.

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning
	Read This Manual/Consult Instructions Before Use
	Caution/consult accompanying documents before use
	Use by Date
	Date of Manufacture
	Batch Code
	Catalogue Number
	Serial Number

Symbol	Meaning
	Manufacturer
	Temperature Limitation
	Humidity Limitation
	Non-sterile
	Fragile, Handle with Care
	Keep Device Dry
	Product Packaging is Recyclable

ELECTROMAGNETIC COMPATIBILITY

Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Emmission Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	This monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

ELECTROMAGNETIC COMPATIBILITY

Table 2

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this monitor requires continued operation during power mains interruptions, it is recommended that this monitor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ELECTROMAGNETIC COMPATIBILITY

Table 3

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING Guidance and manufacturer's declaration - electromagnetic immunity

This monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of this monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{\frac{P}{f}}$ $d = 1.2 \sqrt{\frac{P}{f}} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{\frac{P}{f}} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,** should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this monitor is used exceeds the applicable RF compliance level above, this monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this monitor.

** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC COMPATIBILITY

Table 4

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING
Recommended separation distances between portable and mobile RF communications equipment and this monitor

This monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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PRODUCT SPECIFICATIONS

NAME	Tens Pain Management Solution
MODEL NUMBER	25517 (VH 56-032)
CHANNEL	Single Channel
OUTPUT VOLTAGE	0~34V(vp) +/- 20% (at 500 Ohm load)
PULSE INTENSITY	0~60mA +/- 20%; 0~15 stages Adjustable.
PULSE WIDTH	50~100 uS
PULSE RATE	0 ~ 200Hz
AUTOMATIC SHUT-OFF	After 20 minutes of operation upon completion of therapy treatment
POWER SOURCE	3V CR2032
OPERATION ENVIRONMENT	Temperature 50 F - 104 F (10 C - 40C); Humidity 40% - 90%
STORAGE ENVIRONMENT	Temperature 14 F - 140 F (-10 C - 60 C); Humidity 30% - 95%
DEVICE DIMENSIONS	2-3/4" x 4-1/2"x1/3"
WEIGHT	Controller weight 6.7 oz (with battery)
INCLUDED	Control unit, 1 electrode gel pads set, battery compartment opener, detailed guidebook, replacement electrodes order form, one CR2032 battery
AVAILABLE SEPARATELY	Replacement electrodes 3-pack, Model 25518 (VH 22-037)

Specifications are subject to change without notice

1-YEAR LIMITED WARRANTY

Congratulations on your purchase of the Tens Pain Management Solution. Your Tens Pain Management Solution is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warrantor warrants that its Verde Tens Pain Management Solution will be free from defects in materials and workmanship under normal consumer usage for a period of one year for the original purchaser of the product.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Tens Pain Management System for commercial purposes or subjecting the Tens Pain Management System to abnormal usage or conditions; or (d) other acts which are not the fault of the Warrantor, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Tens Pain Management Solution.

If the Tens Pain Management Solution fails to conform to this limited warranty, return the Tens Pain Management Solution postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any

1-YEAR LIMITED WARRANTY

damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warrantor will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTOR BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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