caring mill™

# Instruction Manual

Deluxe Automatic Wrist  
Blood Pressure Monitor

Model BP3NC1-2AFSA

## Important Product and Safety Information

---

**Follow Instructions for Use.** This document provides important product operation and safety information regarding this **Blood Pressure Monitor**. Please read this document thoroughly before using the device and keep for future reference.

This blood pressure monitor is an automatic digital blood pressure measuring device for use by adults on the wrist at home or in the doctor's office. It enables a very fast and reliable measurement of the systolic and diastolic blood pressure, as well as the pulse rate, by way of the oscillometric method. This device detects the appearance of irregular heartbeats during measurement and provides a warning signal when the irregular heartbeat is detected.

### **This blood pressure monitor IS intended to be used:**

- For self-measurement/monitoring of blood pressure and pulse in adults; whereas the person being measured may be the user/operator of the device.
- Within a home healthcare environment.
- With a cuff located upon the user's wrist; ½ inch from the hand.

### **This blood pressure monitor IS NOT intended for use with:**

- Pregnancy
- Preeclampsia
- Children under 12
- Neonatal patients

**WARNING – Potentially hazardous situation that if not avoided may result in serious injury or death.**

- Self-measuring means monitoring, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you independently alter the dosages of any drugs prescribed by your physician.

- Consult your physician before using this device if any of the following or similar conditions are present: arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, preeclampsia, renal diseases. Motion during measurement, including trembling or shivering may affect the measurement.
- Though not for use with children under 12, ensure that any children around this device are supervised; some parts are small enough to be swallowed.

### **DO NOT:**

- × Use this device if you think it is damaged or anything appears unusual (e.g. sporadic operation, open packaging upon purchase).
- × Use the displayed pulse for checking the frequency of heart pacemakers as this device is not suitable for this action.
- × Open/modify this device; inaccuracy and/or harmful injury may result.
- × Conduct many frequent measurements as this may result in harmful injury due to blood flow interference.
- × Place the Cuff over a wound as this may cause further injury.
- × Place and pressurize the Cuff over/near any present intravascular access or therapy, or arteriovenous shunt, as this may cause blood flow interference and result in harmful injury.
- × Place and pressurize the Cuff over a limb near the side of a mastectomy as this may cause harmful injury.
- × Use the Blood Pressure Monitor on a limb simultaneously with other medical equipment on the same limb due to possible interference with such medical equipment.
- × Maintain pressure in the Cuff applied to the limb for a prolonged amount of time. Ensure that circulation in the limb is not impaired by checking circulation if prolonged/repeated exposure to pressure occurs.

- × Service the device when being used or when power is supplied. When the device is not in use and power is removed, batteries may be replaced by the user. No other parts/components are accessible.

**CAUTION – Potentially hazardous situation that if not avoided may result in minor/moderate injury, property damage, and/or damage to the device**

- This device contains sensitive electronics components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy. Move the device to another location if interference is determined.
- When not using the Blood Pressure Monitor for extended lengths of time, remove the batteries to avoid potential battery leakage and damage to the monitor.
- When replacing the batteries, ensure all the batteries are replaced simultaneously to avoid battery damage and potential damage to the monitor. Microlife USA does not recommend using rechargeable batteries.
- To avoid inaccurate measurements and to lessen any discomfort from Cuff pressure, ensure the Cuff is placed correctly on the limb and fits correctly when snug (not tight), as indicated by markings with the Cuff.
- Consult your physician in cases of frequent irregular heartbeat detections.
- This Risk Indicator feature is provided in order to help you understand your potential blood pressure risk. However, this feature is neither a diagnosis nor a substitute for a medical examination. It is important to consult with your physician to determine your risk.

**DO NOT:**

- × Drop this device or expose it to strong vibrations; sensitive components may be affected resulting in inaccuracies and/or operational issues.

- ✗ Use the Blood Pressure Monitor outside of its specified operation temperature and humidity rating, or if stored outside of its specified storage temperature and humidity rating. Avoid storage in direct sunlight.
- ✗ Use this device in a moving vehicle; inaccurate measurements may result.
- ✗ Use third party accessories. Only use Microlife authorized accessories, such as cuffs or AC adapters, as those not approved for use with the device may provide inaccurate measurements, injury, and/or damage the device.

## STANDARDS

This medical device is compliant with:

- Medical device and non-invasive blood pressure monitor standards IEC 60601-1, IEC 60601-1-2, IEC 60601-1-11, and AAMI/ANSI/IEC 80601-2-30, and
- Electromagnetic standards IEC 60601—1-2, along with FCC Part 15, and
- Clinical Testing per standard ISO 81060-2:2013 was conducted on blood pressure device using the same measurement technology.

Please note: According to international standards, your monitor should be checked for accuracy every 2 years.

## TYPE

IP20: Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Keep Dry. 

Type BF Applied Part

## Expected Life

Monitor 5 Years  
Cuff 2 Years

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

## **FCC**

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications to the product are not approved by Microlife USA and could void the user's authority to operate the equipment under FCC jurisdiction.

This equipment has been verified to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by increasing the distance between the product and the affected device; or if applicable, reorientation the receiving antenna, or moving the product's power plug to another receptacle.

# Wrist Blood Pressure Monitor

## Instruction Manual

### Table of contents

---

#### 1. Introduction

---

- 1.1. Features of the Wrist Blood Pressure Monitor
- 1.2. Gentle+ Measurement System
- 1.3. Important information about self-measurement

#### 2. Important information about blood pressure and its measurement

---

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?

#### 3. Components of your blood pressure monitor

---

#### 4. Using your blood pressure monitor for the first time

---

- 4.1. Installing the batteries
- 4.2. Setting the date and time
- 4.3. Select the User

#### 5. Taking a blood pressure measurement

---

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Discontinuing a measurement
- 5.6. Irregular Heartbeat Detection
- 5.7. Traffic Light Indicator
- 5.8. Memory - Displaying the last 99 measurements

#### 6. Battery Charge Indicator

---

**7. Error messages/malfunctions**

---

**8. Care and maintenance**

---

**9. Limited Warranty**

---

**10. Technical specifications**

---

**11. How to contact us**

---

## **1. Introduction**

---

### **1.1. Features of the Wrist Blood Pressure Monitor**

This blood pressure monitor with Hypertension Risk Indicator is a fully automatic, digital blood pressure measuring device for use on the wrist. This monitor enables very fast and reliable measurement of your systolic and diastolic blood pressure as well as pulse by way of the oscillometric method of measuring. This device offers a very high and clinically tested measurement accuracy and is designed for maximum user-friendliness. Before using it, please read through this instruction manual carefully and keep it in a safe place. For additional questions on the subject of blood pressure and its measurement, please contact your doctor.

### **1.2. Gentle+ Measurement System**

Gentle+, our most advanced system to date, combines clinical accuracy and premium features including:

- A quieter, more comfortable measurement experience that reduces the amount of time your arm is under pressure. Less pressure = more comfort!

## **Attention!**

### **1.3. Important information about self-measurement**

Self-measurement is a means of monitoring blood pressure. It should not be confused with diagnosis or treatment. **You should never alter the dosages of any medication without direction from your physician.** The pulse display is not suitable for checking the functioning of heart pacemakers!!

**In cases of irregular heartbeat, measurements made with this instrument should only be evaluated after consultation with your doctor.**

Electromagnetic Interference: This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens, etc.). These can lead to temporary impairment of the monitor's measuring accuracy.

## **2. Important information about blood pressure and its measurement**

### **2.1 How does high/low blood pressure arise?**

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is affected by way of fine muscles in the blood vessel walls. Your level of arterial blood pressure changes periodically during heart activity. During the "blood ejection" (Systole), the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value). Your blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

## 2.2 Which values are normal?

Blood pressure is too high if your diastolic pressure is above 100 mmHg and/or your systolic blood pressure is over 160 mmHg, while at rest. In this case, please consult your physician immediately. Long-term values at this level endanger your health through continual damage to the blood vessels in you body. Should your systolic blood pressure values lie between 140 mmHg and 159mmHg and/or your diastolic blood pressure values lie between 90 mmHg and 99mmHg. Consult your physician.

Regular self-checks will be necessary. If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician. Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately. If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. Never use the results of your measurements to independently alter the drug doses prescribed by your physician. Which values are normal?

The following standards for assessing high blood pressure in adults have been established by the U.S. National Institutes of Health JNC7, 2003.

<b>Category</b>	<b>Systolic (mmHg)</b>	<b>Diastolic (mmHg)</b>
<b>Normal</b> Self-Check	<120	and <80
<b>Pre-Hypertension</b>	120 - 139	or 80 - 89
<b>Stage 1</b> Hypertension	140 - 159	or 90 - 99
<b>Stage 2</b> Hypertension	≥160	or ≥100

## **Additional information**

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension." Consult your physician.
  - Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.
- \* This device is not intended for use on children 12 years of age or younger. Please ensure that children are supervised around this device, some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.**

## **2.3 What can be done if regular high/low values are obtained?**

- a) Consult your physician.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constrictions caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. Secondary hypertension can be ascribed to specific organ malfunctions. Please consult your physician for information about the possible origins of your own increased blood pressure values.
- d) There are measures you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

### **1) Eating habits:**

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

### **2) Previous illnesses**

- Consistently follow all medical instructions for treating illness such as: diabetes (diabetes mellitus or sugar diabetes), fat metabolism disorder, and gout.

### **3) Habits**


- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea and chocolate).

#### **4) Your physical condition**

- After a medical examination, and with your doctor's approval and direction; exercise.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your physician before beginning your exercise routine. You must receive advice regarding the type and extent of exercise that is appropriate for you.

### 3. Components of your blood pressure monitor

The illustration shows your blood pressure monitor, consisting of:

Fits wrist size 13.5-19.5cm (5.3"-7.7")

## 4. Using your Blood Pressure Monitor for the First Time

---

### 4.1. Installing the batteries

After unpacking the device, insert the batteries. The battery compartment is located on the left side of the device.


- a) Slide the battery cover open.
- b) Insert 2 new batteries (size AAA 1.5V) following the indicated polarity.
- c) If the battery warning appears in the display the batteries are discharged and must be replaced.

#### **Attention!**

- After the low battery indication warning appears, the device will not function until the batteries have been replaced.
- Please use “AAA” Long Life or Alkaline 1.5V batteries. When replacing batteries, replace all at the same time, and ensure they are correctly connected, as shown on the symbols in the compartment.
- We do not recommend using rechargeable batteries with this device.
- If the blood pressure monitor is left unused for long periods, remove the batteries from the device.

## 4.2. Setting the date and time


1. After the new batteries are activated, the year number flashes in the display. You can advance the year by pressing the "M" button. To confirm and then set the month, press the time button.
2. You can now set the month using the "M" button. To confirm and then set the day, press the time button.
3. Please follow the instructions above to set the day, hour and minutes.
4. Once you have set the last minute and pressed the time button, the date and time are set and the time is displayed.
5. If you want to change the date and time, press and hold the time button down for approx. 3 seconds until the user icon starts to flash. Release and repress the button until the year icon blinks. Now you can enter the new values as described above.

Memory Button

Time/Date Button

## 4.3. Select the User

This blood pressure monitor is designed to store 99 measurements for each of two users. Before taking a measurement, be certain that the correct user has been selected.

- a) With the unit off, press and hold the TIME button for 3 seconds until the user icon in the upper left corner of the LCD screen flashes.
- b) Press the "M" button to toggle between users.
- c) Press the START button to make your selection.

## 5. Taking a Blood Pressure Measurement

---

### 5.1. Before the measurement:

- Avoid eating, smoking as well as all forms of exertion directly before your measurement. These factors influence the measurement result. Relax by sitting in an armchair in a quiet atmosphere for about 5 minutes before your measurement.
- Always measure on the same wrist (normally left).
- Compare readings at the same time of day, since your blood pressure changes during the course of the day.

### 5.2. Common sources of error:


#### Note:

Comparable blood pressure measurements always require the same conditions. These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood pressure will be measured! (Each 15cm/6" difference in height results in a measurement error of 10mmHg!)
- A loose cuff causes false measurement values.
- With repeated measurements, blood can accumulate in the respective arm, which can lead to false results. Correctly executed blood pressure measurements should therefore first be repeated after a 1 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away.

### 5.3. Fitting the cuff

a) Remove all objects and jewelry (e.g. wristwatch) from the wrist in use. This is normally the left wrist. Put the cuff over the wrist with the monitor face on the inside of the wrist.

- b) The distance between the cuff and the hand should be approx. 10 mm (1/2").
- c) Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight. No space should remain between the cuff and the wrist.
- d) Lay your arm on a table, with the palm upwards. Support the arm a little with a cushion, so that the cuff rests at about the same height as the heart.


e) Remain seated quietly for five minutes before you begin the measurement.

## 5.4. Measuring procedure

After the cuff has been appropriately positioned the measurement can begin. Remain still during your measurement, do not flex muscles.

- Remain seated quietly for 5 minutes before you begin the measurement.
- Press the ON/OFF button to begin. After pressing the ON/OFF button the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- When the device has detected the pulse, the heart symbol in the display begins to blink.
- After reaching the inflation pressure, the pump stops and the pressure will release. The cuff pressure (large characters) is displayed during the measurement.
- When the measurement has been concluded, the measured systolic and diastolic blood pressure values and pulse will be displayed. (Pulse = number of heart-beats each minute) Example (Fig.): Systolic 118, Diastolic 73, Pulse 75. The measurement results are displayed until you switch the device off or until the device automatically switches off.

## 5.5. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the ON/OFF button can be pressed at any time. The device then immediately lowers the cuff pressure automatically and enters into sleep mode.

## 5.6. Irregular Heartbeat Detection

This symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily), we advise you to tell your doctor. Please show your doctor the following explanation:

### **Information for the doctor on frequent appearance of the Irregular Heartbeat Symbol**

This instrument is an oscillometric blood pressure monitor device that also analyzes pulse frequency during measurement. The instrument is clinically tested.


If pulse irregularities occur during measurement, the irregular heartbeat symbol is displayed after the measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

## 5.7. Traffic Light Indicator


The bars on the left-hand edge of the display show you the range within which the indicated blood pressure value lies. Depending on the height of the bar, the readout value is either within the normal (green), borderline (yellow) or danger (red) range.

The classification is based on standards established by the National Institutes of Health JNC7, 2003. Refer to the chart in section 2.2 of this instruction manual for details of the classifications.

The traffic light bar raises according to your measurement.

- If your measurement has one bar, your measurement is in the green zone, or "Normal" according to NIH standards.
- If your measurement has 2-3 bars, it is in the yellow zone, or "Pre-Hypertension" according to NIH standards.
- If your measurement has four bars, it is in the orange "Stage 1 Hypertension" zone.
- If your measurement has 5-6 bars, it is in the red "Stage 2 Hypertension" zone.

## **5.8. Memory – Displaying the last 99 Measurements**

At the end of a measurement, this monitor automatically stores each result with date and time. This unit stores 99 memories for each of two users.

### **Viewing the stored values**

With the unit off, press the "M" button. The display first shows "A", then shows an average of all measurements stored in the unit.

Pressing the "M" button again displays the previous value. To view a particular stored memory, press and hold the "M" button to scroll to that stored reading.

### **Memory full**

When the memory has stored 99 results for the selected user, a new measured value is stored by overwriting the oldest value.

### **Clear all values**

If you are sure that you want to permanently remove all stored values, hold down the "M" button (the instrument must have been switched off beforehand) until "CL" appears and then release the button. If you do not want to clear the values, press the ON/OFF button. To permanently clear the memory, press the "M" button while "CL" is flashing.

Individual values cannot be cleared.

## 6. Battery Charge Indicator

---

### Batteries almost discharged

When the batteries are approximately 75% used the battery symbol will flash as soon as the instrument is switched on (if at least one of the batteries still has some charge). Although the instrument will continue to measure reliably, you should obtain replacement batteries.

### Batteries discharged– replacements required

When the batteries are discharged, the battery symbol will appear unblinking as soon as the instrument is switched on. You cannot take any additional measurements and must replace the batteries.

1. Open the battery compartment by pushing inward at the two arrows and pulling it out.
2. Replace all of the batteries at the same time. Ensure they are correctly connected, as shown on the symbols in the compartment.
3. The memory retains all stored measurements, however the date, time and the alarms must be reset. For ease of use once the batteries are replaced, the year setting will flash automatically.
4. To set date and time, follow the procedure described in Section 4.2.

### Which batteries and which procedure?

Use two new, longlife 1.5V AAA batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period the batteries should be removed.

**Rechargeable batteries are not recommended for use with this device.**

## 7. Error messages/Troubleshooting

---

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

<b>Error No.</b>	<b>Possible cause(s)</b>
<b>ERR 1</b>	No pulse was detected. Ensure the cuff is properly placed on the inside of your wrist and repeat a measurement
<b>ERR 2</b>	Unnatural pressure impulses influenced the measurement result. Reason: The arm was moved during the Measurement (Artefact). Repeat, keeping still and quiet.
<b>ERR 3</b>	Inflation of the cuff takes too long. The cuff is not correctly seated. Re-position cuff and repeat
<b>ERR 5</b>	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
<b>HI</b>	The cuff pressure is too high. Relax for 5 minutes and repeat the measurement.*
<b>LO</b>	The pulse is too low (less than 40). Repeat the measurement.*

---

**\* If this or any other problem occurs repeatedly, please consult your doctor.**

## Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

### Malfunction

### Remedy

---

The display remains blank when the instrument is switched on.

1. Check batteries for correct polarity.
2. If the display is unusual, re-insert batteries or exchange them.

---

The device frequently fails to measure blood pressure values, or the values measured are too low (too high).

1. Check the positioning of the cuff.
2. Measure blood pressure again in peace and quiet carefully following the details in Section 5.

---

Every measurement produces a different value.

1. Please read the following information listed in Section 5.2  
**"Common sources of error."**  
Repeat the measurement.  
**Please note:** Blood pressure fluctuates continually so successive measurements will show some variability.

---


Blood pressure values measured at your doctors office differ from your measurements.

1. Record the daily development of the measured values and consult your doctor.  
**Please note: Individuals visiting their doctor frequently experience anxiety which can result in a higher reading than at home under resting conditions.**

## 8. Care and maintenance

---

- a) Do not store the unit under direct sunlight, at a high temperature, or in high humidity or dust. Performance may be degraded.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!
- d) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- e) Never open the device! Otherwise the manufacturers warranty becomes invalid.

## 9. Limited Warranty

---

Your Automatic Blood Pressure Monitor is warranted **5 years** by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase.

The 5 year warranty applies to the monitor and the following accessories: cuff. Batteries are not covered by this warranty. There are no user serviceable parts inside.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents.

Professional use, not following the operating instructions, and alterations made to the monitor or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A monitor or accessory as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A monitor or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost.

Please use the below customer service contact information to reach Microlife USA Inc. regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

## 10. Technical specifications:

---

Weight:	148 g (with batteries + cuff)
Size:	85 x 77 x 75 mm (including cuff)
Storage temperature:	-20 to +55°C (-4°F to +131°F)
Humidity:	15 to 90% relative humidity noncondensing
Operation temperature:	10 to 40°C (50°F - 104°F)
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS:	60~255 mmHg
DIA:	40~200 mmHg
Pulse:	40 to 199 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Automatically stores the last 99 .....measurements for two users (total 198)
Measuring resolution:	1 mmHg
Accuracy:	Pressure within $\pm 3$ mmHg or 2% of reading .....>200mmHg, Pulse $\pm 5$ % of the reading
Power source:	2 batteries size AAA 1.5V
Accessories:	storage case
Technical alterations reserved!	
Made in China	

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